GELITA TUFT-IT®

Absorbable Non-Woven Gelatin Hemostat

• Instructions for use

CE 0344

- (EN) Absorbable Non-Woven Gelatin Hemostat p. 04
- (AR) شاش جيلاتيني لانسيجي ماص للجاء p. 07
- (FR) Hémostatique de gélatine non tissé résorbable p. 10
- (DE) Resorbierbare Vlieseinlage auf Gelatinebasis p. 14
- (IT) Emostatico in gelatina assorbibile non in tessuto p. 17
- (PT) Hemostático de gelatina absorvível não tido p. 21
- (RU) Рассасывающееся нетканое
кровоостанавливающее средство из желатина p. 25
- (ES) Hemostático de gelatina absorbible Fibrilar p. 29
- (TR) Emilebilir Nonwoven Jelatin Kanama Durdurucu p. 32
- (NL) Absorbeerbare, non-woven, gelatine hemostaat p. 36

GELITA®
MEDICAL

Complete Care! Completely Safe!

(EN) Explanation of symbols (AR) مدلول الرموز (FR) Signification des symboles
 (DE) Symbolerklärung (IT) Spiegazione dei simboli (PT) Explicação dos símbolos
 (RU) Объяснение символов (ES) Explicación de los símbolos
 (TR) Sembollerin açıklaması (NL) Verklaring van de symbolen

(EN) Fibrillar (AR) لِيُفِّي (FR) Fibrillaire (DE) Fibrillar (IT) Fibrillare (PT) Fibrilar (RU) Состоит из волокон (ES) Fibrilar (TR) Fibriler (NL) Fibrillerend

(EN) Manufacturer (AR) الشركة المصنعة (FR) Fabricant (DE) Hersteller (IT) Produttore
 (PT) Fabricante (RU) Производитель (ES) Fabricante (TR) Üretici (NL) Fabrikant

(EN) Date of manufacture (AR) تاريخ التصنيع (FR) Date de fabrication (DE) Herstellungsdatum
 (IT) Data di produzione (PT) Data de fabrico (RU) Дата изготовления (ES) Fecha de fabricación
 (TR) Üretim tarihi (NL) Productiedatum

(EN) Expiry Date (AR) تاريخ انتهاء الصلاحية (FR) Date de péremption (DE) Verfallsdatum
 (IT) Data di scadenza (PT) Data de validade (RU) Дата истечения (ES) Fecha de caducidad
 (TR) Son Kullanma Tarihi (NL) Vervaldatum

LOT

(EN) Batch code (AR) كود التشغيل (FR) Code du lot (DE) Chargennummer (IT) Codice lotto
 (PT) Código do lote (RU) Код пакета (ES) Número de lote (TR) Seri kodu (NL) Batchcode

REF

(EN) Order number (AR) رقم الطلب (FR) Numéro de commande (DE) Bestellnummer
 (IT) Numero d'ordine (PT) Número de encomenda (RU) Номер для заказа (ES) Número de pedido
 (TR) Sipariş numarası (NL) Ordernummer

(EN) Store dry (AR) تخزين جاف (FR) Conserver au sec (DE) Trocken aufbewahren
 (IT) Conservare in luogo asciutto (PT) Armazenar em local seco (RU) Хранить в сухом месте
 (ES) Almacenar en seco (TR) Kuru saklayın (NL) Droog bewaren

(EN) Keep away from sunlight (AR) لا تعرضه لأشعة الشمس (FR) Ne pas exposer aux rayons du soleil
 (DE) Von Sonnenlicht fernhalten (IT) Conservare al riparo dalla luce solare (PT) Manter protegido da luz solar
 (RU) Хранить вдали от солнечного света (ES) Mantener lejos de la luz solar
 (TR) Güneş ışığından uzak tutun (NL) Beschut tegen het zonlicht bewaren

 5°C **EN** Temperature limitation **AR** حد درجة الحرارة **FR** Limites de température
DE Temperaturbegrenzung **IT** Limite di temperatura **PT** Limitação de temperatura
RU Ограничение температуры **ES** Límites de temperatura **TR** Muhafaza edilmesi gereken sıcaklık aralığı **NL** Temperatuurbeperking

 **EN** Single use **AR** الاستخدام مرة واحدة **FR** Usage unique **DE** Nicht wiederverwenden
IT Monouso **PT** Descartável **RU** Одноразовое использование **ES** Un solo uso
TR Tek kullanımlık **NL** Eenmalig gebruik

 **EN** Do not use if the packaging is damaged **AR** لا يُستخدم في حالة تلف التغليف
FR Ne pas utiliser si l'emballage est endommagé **DE** Bei beschädigter Verpackung nicht verwenden
IT Non utilizzare se la confezione è danneggiata **PT** Não utilizar se a embalagem estiver danificada
RU Не использовать, если упаковка повреждена **ES** No usar si el embalaje está dañado
TR Ambalaj hasarlı ise kullanmayın **NL** Niet gebruiken als de verpakking is beschadigd

STERILE **EN** Radiation sterilization **AR** تعقيم بالإشعاع **FR** Stérilisation par irradiation
DE Sterilisation durch Bestrahlung **IT** Sterilizzazione con radiazione **PT** Esterilização por radiação
RU Стерилизация облучением **ES** Esterilización mediante radiación
TR Radyasyonla sterilizasyon **NL** Stralingssterilisatie

 **EN** Do not re-sterilize **AR** لا تعيد تعقيمه **FR** Ne pas restériliser **DE** Nicht erneut sterilisieren
IT Non risterilizzare **PT** Não esterilizar novamente **RU** Не стерилизовать повторно
ES No reesterilizar **TR** Yeniden sterilize etmeyin **NL** Niet opnieuw steriliseren

 **EN** Use dry for best performance **AR** استخدم جافاً للحصول على أفضل أداء
FR Utiliser à sec pour une performance optimale **DE** Trocken anwenden für optimale Leistung
IT Utilizzare asciutto per le migliori prestazioni. **PT** Utilizar o produto seco para um melhor desempenho
RU Для лучшего результата использовать сухим **ES** Utilizar en seco para un mayor rendimiento
TR En iyi performans için kuru kullanın **NL** Droog gebruiken voor de beste prestaties

 **EN** Consult instructions for use **AR** الرجوع لتعليمات الاستخدام **FR** Consulter les instructions d'utilisation
DE Gebrauchsanweisung beachten **IT** Consultare le istruzioni per l'uso
PT Consultar instruções de utilização **RU** См. инструкцию по применению
ES Consultar las instrucciones de uso **TR** Kullanım talimatlarına başvurun
NL De gebruiksaanwijzingen raadplegen

GELITA-TUFT-IT®

EN

Absorbable Non-Woven Gelatin Hemostat

1. Product description

GELITA TUFT-IT® is a sterile, water insoluble, biodegradable multi-layer non-woven fabric, absorbable hemostat made of 100 % purified porcine gelatin. It is manufactured out of gelatin fibres in a special spinning process. The fleece is intended for use as a topical haemostat as an adjunct to haemostasis by tamponade effect. Due to the high flexibility of the material it can be used to control bleeding over a large area, in surface applications, irregularly shaped bleeding sites or in areas that are difficult to access; the material can be tufted or peeled in layers as required. The fibers continue to adhere to one another; there is no uncontrolled dispersal over the wound site. In use GELITA TUFT-IT® swells minimally as blood and other fluids are absorbed.

Every hemostat is packaged individually in double blister, sterilized by means of gamma irradiation, and is for single use only.

2. Indications

GELITA TUFT-IT® absorbable hemostat for use as an adjunct to hemostasis by tamponade effect, in particular where control of capillary and venous bleeding by pressure, ligature and other conventional procedures is either ineffective or impractical.

3. Contraindications

- Do not use GELITA TUFT-IT® in patients with known allergies/hypersensitivity to porcine products.
- GELITA TUFT-IT® should not be used in closure of skin incisions, as it may interfere with healing of the skin edges. This is due to the mechanical interposition of the gelatin and not to intrinsic interference with wound healing.
- The product should not be used in infected wounds without antibiotics.
- GELITA TUFT-IT® should not be used in conjunction with methyl methacrylate adhesives, for example in orthopaedic surgery, because their presence may reduce the adhesive strength of the bonding agent to the bone.
- GELITA TUFT-IT® should not be used for controlling postpartum bleeding or menorrhagia.
- GELITA TUFT-IT® should not be used in bleeding from large arteries, as this may result in embolization of the blood vessel.

4. Side effects

There have been no reported adverse reactions for this kind of product when used correctly according to these Instructions for Use and when overpacking is avoided.

The formation of tissue granulation following middle ear interventions was observed for this kind of product in at least one animal study.

5. Directions for use

Sterile technique should be applied when using GELITA TUFT-IT®.

GELITA TUFT-IT® should be applied dry and should not be moistened prior to application. The contact with saline or antibiotic solutions before hemostasis is achieved will cause a reduced hemostatic performance. If used in combination with liquid agents, fluids must be added after application of GELITA TUFT-IT® to the wound surface.

Dry gloves and surgical instruments (forceps, scissors) should be used to handle, cut and apply GELITA TUFT-IT®.

GELITA TUFT-IT® is best employed by loose application against the bleeding surface, where it should adhere. Press the product carefully and gently to the wound surface until hemostasis is achieved.

GELITA TUFT-IT® can be cut to fit for use in endoscopic procedures.

The material can remain at the site once the bleeding has stopped. To facilitate bio-degradation and minimize the possibility of foreign body reaction and encapsulation of fluid, only use minimal amounts of the device necessary for hemostasis and remove any excess of material prior to surgical closure. When used correctly and in small amounts, GELITA TUFT-IT® is fully absorbed at the implantation site. Depending on the application procedure and the area of surgery, GELITA TUFT-IT® undergoes biodegradation within four weeks. However, bio-degradation depends upon several factors including the amount used, degree of saturation with

blood and the tissue bed. Any excess material not removed at the time of surgery may present as a recurring mass or lesion or postoperative abscess formation upon imaging.

Any excess material must be discarded in the hazardous dust bin of the Operation Room.

Important notice!

GELITA TUFT-IT® expands minimally with the absorption of fluids. In radical cavities; laminectomy procedures; around or in the proximity to foramina in bone; areas of bony confine; the spinal cord; the optic nerve and chiasm; or closed tissue spaces with presence of bone, GELITA TUFT-IT® should be removed when bleeding has stopped, if possible. Failure to do so might result in unintended pressure on neighbouring structures which may cause pain for the patient or could create the potential for nerve damage.

Safety of single or multiple applications of the GELITA TUFT-IT® at the (same) bleeding area has not been investigated in humans.

While most reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures.

By absorbing fluid GELITA TUFT-IT® may expand and impinge on neighbouring structures. Therefore, care should be taken to avoid over-packing, particularly within bony cavities, as it may lead to interference with normal function and may eventually lead to compression necrosis of surrounding tissue.

In cases of a postoperative infection, re-op-

eration may be necessary to remove infected material and allow drainage.

General warnings

GELITA TUFT-IT[®] is not intended as a substitute for careful surgery and the proper use of sutures and ligatures or other conventional procedures for hemostasis.

Avoid using the product in infected areas. Embedding the product in a contaminated wound without drainage may lead to complications and should be avoided.

There is no own clinical data available for the application of the product in children and during pregnancy. Therefore GELITA TUFT-IT[®] should only be used in those cases if the benefit of its application clearly exceeds the risks.

GELITA TUFT-IT[®] is not recommended for the primary treatment of coagulation disorders.

For application in patients with coagulation disorders no scientific data is available.

The use of the product in patients on anticoagulant therapy may cause an extended time to hemostasis.

Care should be taken to avoid the introduction of fragments of the material directly into the circulatory system.

GELITA TUFT-IT[®] should not be used for embolization or to treat bleeding from large arteries, as hemostatic devices used near moderate to large blood vessels may result in embolization of the blood vessel. Such embolization has been associated with severe adverse effects, including fever, duodenal and pancreatic infarct, and embolization of lower extremity vessels, pulmo-

nary embolization, splenic abscess, necrosis, asterixis, and death.

Precautions should be taken in otorhinolaryngology surgery to assure that none of the material is aspirated by the patient. (Examples: controlling haemorrhage after tonsillectomy and controlling epistaxis)

Due to the neutral properties of the product, it can be combined with drugs (for information please read the package insert of the respective medication).

The safety and efficacy of GELITA TUFT-IT[®] in combination with those agents such as topical thrombin has not been evaluated in clinical trials. If in the physician's judgement concurrent use of topical thrombin or other agents are medically advisable, the respective instructions for use should be consulted.

Cell Savers

The use of GELITA TUFT-IT[®] with a cell saver is not contraindicated. However, make sure to use proper technique when using GELITA TUFT-IT[®] in conjunction with a cell saver. As with other hemostatic agents, GELITA TUFT-IT[®] should not be aspirated into autologous blood salvage circuits.

6. General storage and handling

GELITA TUFT-IT[®] should be stored dry at 5 to 30°C in its sealed outer protective packaging. For short periods (e.g. transportation) deviations of +/- 10°C can be permitted. Keep away from direct sunlight.

Remove the product from the sterile package

using sterile technique. Prior to opening, the integrity of the sterile barrier should be confirmed. Discard the product if the package is damaged.

Any unused contents must be discarded.

The packaging is subject to and conforms to European and international legislation and other

applicable standards. The packaging protects the product from external influences and guarantees its sterility during storage and transport, under normal use.

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GELITA MEDICAL GmbH